

FEATHER Skin Grafting Blade

SKU: MP23566



Categories: [Surgical Holloware](#)

Description

We are among the prestigious organizations engaged in importing quality range of Skin Graft Blade from Feather (Japan). These knives are manufactured by highly professionals, using latest machines and high grade material. SPECIFICATIONS:

- Manufactured by using high quality stainless steel
- Handle and stem equipped with mounted and adjustable roller
- Adjustable thickness of graft and blade holder
- Designed to fit mostly all handles
- Available in both carbon and stainless steel supplied sterile
- Packed in a medical grade pouch

USES:

- Used for taking split skin grafts
- Used for coverage of raw surface composed of healthy vascularized or clean granulations to obtain healing and limit deformity and disability

PROCEDURE:

- The limb is held firmly after preparing the donor site
- The knife is held pressed firmly against the skin and the graft is removed with a steady to and for sawing motion
- The desired thickness of the graft can be obtained by adjusting the depth of the knife
- A very thin graft is translucent and the knife blade underneath appears bluish gray in color
- It leaves behind numerous small closely packed or confluent bleeding points
- The subcutaneous fat is seen in case of a full thickness graft