

# Accuflow Drip Monitor

SKU: MP49900



Categories: [Diagnostic Equipments](#), [Health Innovations](#), [Hospital use innovations](#)

## Short Description:

The conventional gravity infusion method requires the nursing staff / caregivers to manually calculate the required flow rate and monitor the infusion at periodic intervals. Given the work load faced by nursing staff in a typical hospital scenario, periodic patient monitoring for IV is an additional burden and often times, patients or their care takers complain of IV not being infused and patient discomfort. Having an IV monitoring makes nurses job easy and efficient and eliminates scope for any human errors or late intervention in case of any alarms.

## Description

The conventional gravity infusion method requires the nursing staff / caregivers to manually calculate the required flow rate and monitor the infusion at periodic intervals. Given the work load faced by nursing staff in a typical hospital scenario, periodic patient monitoring for IV is an additional burden and often times, patients or their care takers complain of IV not being infused and patient discomfort. Having an IV monitoring makes nurses job easy and efficient and eliminates scope for any human errors or late intervention in case of any alarms. Accuflow Drip Monitor is designed to help you in IV monitoring Accuflow Drip Monitor (ADM) device is used to monitor the gravity infusion and to provide alarms during the critical conditions.

### Product Key Features:

- Monitor IV closely
- Precise Flow Measurement
- Critical Alarm display
- Reliable Drop Sensing
- Highly Affordable
- Zero Learning Curve
- Improved Patient safety
- Reduce Nurse Workload
- Remote Monitoring Station with alarm dashboards